



V6 Full Digital Ultrasound Diagnostic System



Dustproof, waterproof, adaptable to harsh environments



Full-viewing angle IPS high-resolution screen for multi-angle viewing




English touch buttons



High-capacity lithium battery, can work continuously for 8 hours

 **SonoStar**



Detection content

- Detection of pregnancy, number of fetuses, foetal development
- Obstetric disease detection, detecting uterine diseases, changes in ovarian follicles and corpus luteum, cysts, pyometra, fibroids and other diseases
- Detect the gestational age and expected date of delivery, distance, circumference, area, heart rate of cattle, horses, donkeys, sheep, pigs, dogs, cats, camels and other animals
- One-click measurement of back fat, eye muscle area, lean meat percentage, etc.

Main performance

- The whole machine is fully waterproof and can be washed and cleaned with water
- Optional video glasses provide excellent viewing experience
- Ultra-high-definition 8-inch large LCD medical Chinese touch display, dual internal and external screens, explosion-proof and scratch-proof
- Full-viewing IPS high-resolution screen, capable of observing the measured target from multiple angles
- 10000 mAh battery, can work continuously for 8 hours
- One-click smart screen rotation, images and touch buttons can be flipped, making it easy for both left and right hands to operate
- Using metal plate radiator, excellent heat dissipation function
- Automatic backfat measurement function (intelligent one-button measurement)
- Supports a variety of electronic probes. Each probe has five-stage frequency conversion function. The probe does not have a conversion box, making it more convenient to use.
- Multiple movie playback functions, any image measurement annotation can be selected and saved
- 32g internal storage space, massive images can be stored at will
- It has grid ruler function and dot ruler function, which can quickly read the size of the measured target.
- A rich variety of animal-specific software measurement packages that can automatically calculate pregnancy cycles, accurately and conveniently
- Sixteen colors are adjustable, making it easy to observe the location of lesions
- Waterproof interface: supports USB interface, HDMI interface, TE interface, power interface
- It has settings such as time and date, video output, TV format, key tone switch, parameter mode, standby protection, multi-language selection, customized gestational age, etc.
- One-click restore factory parameters function, the operator does not need to worry about parameter adjustment chaos

Excellent Images



Calf embryo



Calf trunk



Calf teeth



Calf Bone



Optional video goggles for ultimate viewing experience



Strap-type design makes the operator more convenient and comfortable



Environmentally friendly and lightweight Aluminum-plastic packaging case

Guangzhou Sonostar Technologies Co., Limited

Address: C Building, #27 Yayingshi Road, Science Town, Guangzhou, China

Tel/ Fax: +86-20-32382095

Web: www.sonostarmed.com

E-mail: sales@sonostar.net